IVD

Lomina Vitamin D Test, Semi-quantitative
Instructions for Use (IFU)

ENGLISH

[PACKING SPECIFICATION]

Variants	1 piece in a package - box/pouch	25 pieces in a package - box	50 pieces in a package - box
REF	La-D-Pro/1B or 1P	La-D-Pro/25B	La-D-Pro/50B
IVD test strip in a plastic cassette [pcs]	1	25	50
Sterile Lancet [pcs]	1	25	50
Stabilization fluid - buffer [pcs]	1	25	50
15 µl dropper [pcs]	1	25	50
Disinfecting pad [pcs]	1	25	50
IFU [pcs]	1	1	1

[INTRODUCTION]

Vitamin D is a group of fat-soluble secosteroids that increase intestinal absorption of calcium, iron, magnesium, phosphate and zinc. In humans, the most important compounds in this group are vitamin D3 and vitamin D2. Vitamin D3 is naturally produced in human skin by exposure to ultraviolet radiation and vitamin D2 is obtained mainly from food. Vitamin D is transported to the liver where it is metabolized to 25-hydroxy vitamin D. The concentration of 25-hydroxy vitamin D (including D2 and D3) in the blood is considered the best indicator of vitamin D status. Vitamin D deficiency is now considered a worldwide problem. Virtually every cell in our body has receptors for vitamin D, which means they all require "sufficient" levels of vitamin D to function adequately. The health risks associated with vitamin D deficiency are much more serious than previously thought. Vitamin deficiency is associated with a variety of serious diseases: osteoporosis, osteomalacia, multiple sclerosis, cardiovascular disease, pregnancy complications, diabetes, depression, stroke, autoimmune disease, influenza, various cancers, infectious diseases, Alzheimer's disease, obesity and higher mortality rates, etc.

[INTENDED USE]

The Lomina Vitamin D Test is a rapid chromatographic immunoassay for the semi-quantitative detection of 25-hydroxyvitamin D (25(OH)D) in the human body. It is designed for detection from a finger prick sample of whole blood. This test provides an indicative diagnostic result and can be used to screen for vitamin D deficiency.

[TEST PRINCIPLE]

The Vitamin D test is based on the principle of competitive binding immunoassay. During testing, the compound migrates through the membrane where an immunological reaction occurs during capillary migration. The membrane is pre-coated with 25-OH-vitamin D antigen in area (C). The 25-OH-vitamin D in the sample is competitive with the 25-OH-vitamin D in the test area. After a limited amount of anti-25-OH-vitamin D antibodies in the conjugate, all conjugated antibodies are saturated with 25-OH-vitamin D in the sample, no coloured line is formed in the test region, indicating a concentration greater than 30 ng/ml. The presence of a visible coloured line in the test area (T) indicates the concentration of 25-OH-vitamin D at a level below 30 ng/ml. The higher the concentration of 25-OH-vitamin D in the sample, the lighter the line will be in this manual. As a check on the procedure, a colored line C will always appear in the control area to indicate that the correct volume of sample has been added and that membrane activation has occurred.

[PACKAGE CONTENTS]

- IVD test strip in a plastic cassette.
- Desiccant.
- Lancet for dispensing blood sample collection (HTL-Strefa S.A. C E0344; Typ Actilance INTROLINE) or wellion MED TRUST C E0197; Type 23G double protection (double-sided).
- Stabilization fluid (buffer).

[CROSS-REACTIVITY]

There was no interference with the following substances: acetoaminophen, acetalsalicylic acid, ascorbic acid, atropine, bilirubin, caffeine, gentisic acid, glucose, hæmoglobin, triglycerides, albumin.

[PERFORMANCE CHARACTERISTICS]

Testing was performed using the PCR Vitamin D Kit. Testing was performed in a European GLP compliant laboratory.

[NOT INCLUDED IN THE PACKAGE]

Stopwatch/clock for measuring the duration of the test.

[ATTENTION!]

- Please read all information in this manual before performing the test.
- For use by medical professional in in vitro diagnostics only.
- Do not eat, drink or smoke in the area where samples or kits are handled.
- Store in a dry place at a temperature of 2-30°C, avoid areas with excessive humidity, the foil packaging is damaged or has been opened, do not use it.
- The test kit is intended for preliminary testing only.
- Strain or observe the specified sample reading time.
- Use the test only once. Do not disassemble or touch the test window of the test cassette.
- The kit must not be frozen or used beyond the expiry date on the packaging.
- Keep out of the reach of children.
- Dispose of the used test in accordance with local regulations.

[STORAGE CONDITIONS]

1. The test pack shall be stored at a temperature of 2°C to 30°C.
2. The components of different batches shall not be mixed. Each component is stable under the specified conditions and may reach the specified shelf life of the kit.

[ACCURACY]

The accuracy of the Lomina Vitamin D test is 93.8%.

[TABLE OF REFERENCE RANGES]

25-OH Vitamin D level	Reference range [ng/mL]	Reference range [nmol/L]
Deficient	0-10	0-25
Insufficient	10-30	25-75
Sufficient	30-100	75-250

[LINE COLOR INTENSITY]

Compare the color intensity of the line in the test area (T) with the reference line below.

[INTERPRETATION OF THE TEST RESULTS]

The accuracy of the Lomina Vitamin D test is 93.8%.

[STABILISATION]

The result is interpreted by the visible line in the "T" area on the body of the test cassette (see figure below). The results of the visual inspection are as follows:

[THE RESULT IS INTERPRETED BY THE VISIBLE LINE IN THE "T" AREA ON THE BODY OF THE TEST CASSETTE (SEE FIGURE BELOW). THE RESULTS OF THE VISUAL INSPECTION ARE AS FOLLOWS:]

[THE ACCURACY OF THE LOMINA VITAMIN D TEST IS 93.8%]

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[LAGERUNGSBEDINGUNGEN]

- Der Test set shall be stored at a temperature of 2°C to 30°C.
- Components of different batches shall not be mixed. Each component is stable under the specified conditions and may reach the specified shelf life of the kit.

[TESTLIMITATIONEN]

1. The Lomina Vitamin D Test provides only a semi-quantitative analytical result. A secondary analytical method must be used to obtain a confirmed result.

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